



# Compatible Hardware and Accessories for PatientTouch® Communications for Android for 2.9 and 4.x

---

Reproduction in any manner whatsoever without the written permission of PatientSafe Solutions, Inc. is strictly prohibited. Changes in equipment, software, or procedures may occur periodically; information describing these changes will be included in future editions of this document. Information in this document is subject to change without notice and does not represent a commitment on the part of PatientSafe Solutions to provide additional services or enhancements.

All patient names and medical histories depicted in these materials are fictitious. Any resemblance to an actual person or case is purely coincidental. All drugs, drug orders, and dosages depicted in screen shots, videos, diagrams, or other media are for illustrative purposes only. PatientSafe Solutions, Inc. makes no recommendation or representation about any treatment or dosage. These materials are subject to revision by PatientSafe Solutions, Inc. from time to time, at any time.

**PatientTouch** is a registered trademark of PatientSafe Solutions, Inc. iPhone® and iPod touch® are registered trademarks of Apple Inc. Other product or company names are the trademarks and/or registered trademarks of their respective owner.

## Device Options and Specifications

The following devices have been certified to be compatible with PatientTouch® Communications for Android 2.9 and above.

Android devices are currently not certified for the PatientTouch Clinical application, however, support for this option is in progress.

All of these devices can be used for shared device and BYOD, however, **PatientSafe Solutions™ is working with vendors to certify specific devices for VoIP.**

Certified Android Platforms: Pie, Oreo, Android 10

PatientSafe Solutions recommends a battery charger and case, however, we have not tested or approved any specific brands.

Device
Google Pixel 2
Google Pixel 3/3XL
Samsung Galaxy S7 Edge
Samsung Galaxy S6
Samsung S9
Samsung S7 Edge
Samsung Galaxy J7
Samsung Galaxy S8+
Samsung S8
Samsung Note 8
Samsung S9
Galaxy S8
LG Nexus 4
Samsung Galaxy S4

Samsung Note 2
Galaxy S9+
Galaxy S7
Galaxy S10
Galaxy S10+

**Note:** The devices in this list are tested and verified by PatientSafe Solutions™. If you are having problems with your PatientTouch® Communications for Android application and you are not using one of the above devices, please switch to a verified device or contact the Technical Support Team:  
Phone: (888) 310-3105 - [support@patientsafesolutions.com](mailto:support@patientsafesolutions.com)

\* 5GHz network support required for VOIP (Voice)