

4.9.2 Server Only Release

This release note is to inform customers of new features, enhancements, and resolved issues in version 4.9.2 of the PatientTouch System.

Release Information

PSS recommends you upgrade your server to version 4.9.2.

Affected Hospital Staff

All PatientTouch 4.x users.

Comms

New Features and Enhancements

Small Group Limit for Video Calls

Description: Default video settings have been updated so users will no longer reach the 4 person small group participant limit. For example, if there were less than four users on a call and they tried to add additional users that surpassed 4 persons, they would not be able to add the additional users without restarting the call. This new feature eliminates the 4 person small group limit for video calls.

Audience: All PatientTouch Clinical and Comms 4.x users.

Impact: Groups of less than 4, are now allowed to add additional users up to and over 4 persons.

Third Party Requirements: N/A

Configuration and Maintenance: N/A

Devices and Applications: PatientTouch Clinical iOS Application, PatientTouch Communications Application (iOS)

Additional Resources: N/A

Global System Release Notes

Resolved Issues

1. Fixed issue where badge count was indicating a certain number of tasks that needed to be completed when in reality there were no tasks to be completed at all.

Specimen Collection Release Notes

Resolved Issues

1. Fixed issue reviving lab orders in Clinical Manager.